Skylar's Run

THYNK
thynk.com

Quick Start Guide for Teachers & Students

Step 1: Charge your Headset and Connect your iPad to Wi-Fi

- Make sure to use a fully charged headset - the light on the headset should turn green!
- Ensure the iPad is connected to your Wi-Fi network.

Step 2: Put on Your Headset Once it is Charged

- Using a piece of blotting paper, blot the skin on your forehead and on the mastoid bone behind your right ear with a piece of blotting paper.
- Unplug the headset from the charger. Place the headset on as modeled during your orientation meeting and adjust as needed. All four sensors should be in contact with your skin.
- If you have any connection issues, add a few drops of saline to each sensor.



Step 3: Play the Game! Train your Brain!

- On your iPad, launch the Skylar's Run app.
- Input your username and password:
  - Username: _____
  - Password: _____
- Select your headset ID when it appears on screen.
- Begin playing Skylar's Run!

Questions?

Email support@thynk.com